

Development and validation of a quantitative confirmatory method for 30 β -lactam antibiotics in bovine muscle using liquid chromatography coupled to tandem mass spectrometry

M. Di Rocco^{a,b}, M. Moloney^a, S. Earley^c, B. Berendsen^d, A. Furey^b, M. Danaher^{a*}

^aFood Safety Department, Teagasc Food Research Centre, Ashtown, Dublin 15, Ireland

^bMass Spectrometry Research Group, Department of Physical Sciences, Cork Institute of Technology, Bishopstown, Cork, Ireland

^cVeterinary Public Health Regulatory Laboratory, Department of Agriculture Food and the Marine, Backweston, Celbridge, Ireland

^dRIKILT Wageningen UR, Department of Veterinary Drug Residues, Wageningen, The Netherlands

Keywords: β -Lactam residue; Penicillin; Cephalosporin; Carbapenem; UHPLC-MS/MS; Bovine muscle

Corresponding author: martin.danaher@teagasc.ie

Abstract

A method was developed for the confirmatory and quantitative analysis of 30 β -lactam antibiotic residues in bovine muscle. The method includes 12 penicillins (amoxicillin, ampicillin, cloxacillin, dicloxacillin, mecillinam, methicillin, nafcillin, oxacillin, penicillin G, penicillin V, piperacillin, ticarcillin), 12 cephalosporins (cefacetrile, cefadroxil, cephalixin, cefalonium, cefazolin, cefoperazone, cefotaxime, cefquinome, cefuroxime, desacetyl cephalirin, desfuroylceftiofur cysteine disulfide, desfuroylceftiofur dimer), five carbapenems (biapenem, doripenem, ertapenem, imipenem, meropenem) and faropenem. Samples were

extracted using a simple solvent extraction with acetonitrile:water (80:20, v/v) and C₁₈ dispersive solid-phase extraction (d-SPE) clean-up, followed by ultra-high performance liquid chromatography coupled to tandem mass spectrometry (UHPLC-MS/MS) detection. Chromatography was performed on a reversed phase CSH C₁₈ column, using a binary gradient separation comprising of 0.01% formic acid and 0.2 mM ammonium acetate in water (mobile phase A) and 0.01% formic acid in acetonitrile (mobile phase B). The mass spectrometer was operated in the positive electrospray ionisation mode (ESI(+)). Validation was performed following the 2002/657/EC guidelines. Trueness ranged between 69% and 143% and precision ranged between 2.0% and 29.9% under within-laboratory reproducibility conditions. Absolute recoveries ranged from 40.3% to 89%. The developed method uses minimal sample preparation and 30 test samples can be analysed by a single analyst in a single day. To the best of our knowledge, this is the first method for carbapenems in foodstuff that does not require derivatisation.

1. Introduction

The β -lactams are key antibiotics used in both human and veterinary medicine [1]. All of these drugs have a common four-membered ring in their molecular structure (Fig. 1). The penicillins, cephalosporins and carbapenems represent the most important β -lactam drug groups. In the penicillin and cephalosporin structures, the four-membered ring is fused to a five-membered thiazolidine or a six-membered dihydrothiazine ring, respectively. The carbapenems are similar to the penicillins, with a sulphur atom replaced by a carbon atom. Faropenem, a carbapenem-related compound, is another important member of the β -lactam class [2].

The β -lactams act by inhibiting the enzymes (PBPs, Penicillin-Binding Proteins) involved in the biosynthesis of the peptidoglycan, a fundamental constituent of the bacterial cell wall [3]. This is particularly important for actively multiplying bacteria, in which the β -lactams can cause cell lysis. Penicillins are classified as broad- or narrow-spectrum antimicrobials, based on their activity, or according to their susceptibility to β -lactamases produced by bacteria [1]. Cephalosporins are less susceptible to the action of β -lactamases and can be categorised into different generations [4]. In particular, their activity against Gram-positive bacteria decreases between the first and third generations, while their effectiveness against Gram-negative bacteria increases. The fourth-generation cephalosporins are broad-spectrum antimicrobials that are active against both Gram-positive and Gram-negative organisms [1].

Penicillins and cephalosporins are often administered parenterally and orally to food-producing animals to prevent or treat bacterial infections. An inappropriate use of antibiotics may lead to residues in food, which can cause health hazards, such as allergic reactions in sensitive individuals. In addition, concerns have been raised over the excessive usage of antibiotics, which is linked to the increase of antimicrobial resistance in livestock and the potential transfer of resistant bacteria and resistant genes to humans and human pathogens, respectively [3]. This is particularly important among the third- and fourth-generation cephalosporins, due to the reliance of these drugs in human medicine as broad-spectrum antibiotics [5, 6]. Consequently, Maximum Residue Limits (MRLs) are established for different species and target tissues [7] and accurate methods become essential to measure the presence of non-compliant antibiotic residues in food.

The monitoring of β -lactams in food of animal origin is complicated because of the differences in authorisation in food producing species. There are nine penicillin and eight

cephalosporin active ingredients listed under Commission Regulation (EU) 37/2010 [7]. Many of the penicillin drugs are licensed for use in multiple species. For example, amoxicillin, ampicillin, cloxacillin, dicloxacillin, oxacillin and penicillin G are allowed in the treatment of all food-producing species. The other penicillins are licensed more specifically, e.g. nafcillin (all ruminant species), penicillin V (poultry and swines) and mecillinam (cattle only).

The usage of cephalosporin antibiotics in food production has attracted much attention in recent years with calls to restrict the applications of third- and fourth-generation drugs. The European Commission recommends that the use of a narrow-spectrum antimicrobial should always be the first choice unless prior susceptibility testing shows that this would be ineffective [8]. Consequently, third- and fourth-generation cephalosporins should only be administered systemically to animals when clinical conditions respond poorly to narrow spectrum antibiotics. Among the first-generation cephalosporins, cefacetrile, cephalixin, cefalonium and cephapirin are approved for the treatment of cattle, while cefazolin is approved for cattle, goats and swines. Cefoperazone and ceftiofur belong to the third generation and are allowed in the treatment of cattle and all mammalian food-producing species, respectively. Cefquinome, a fourth-generation cephalosporin, is approved for cattle, horses and swine [7]. The carbapenems are broad spectrum antibiotics, resistant to most β -lactamases. These drugs are not registered for use in food-producing animals, but play a critically important role in human medicine [9]. As a consequence, analytical methodologies should be capable of detecting MRL breaches for licensed drugs and the presence of residues from unauthorised applications.

Different screening approaches are available for detecting antibiotic residues in food samples (immunoassays, microbial inhibition assays, receptor assays). The main advantages of these tests are simplicity, low cost and short analysis time, with no need for laborious sample preparation. However, most of these methods cannot identify or quantify the antibiotic that is present. In addition, some of these tests can lead to false positive or false negative results [10]. Liquid chromatography coupled to mass spectrometric detection has become the preferred technique to quantify and confirm the presence of β -lactams in foodstuffs. Several LC-MS(/MS) methods have been reported for the analysis of β -lactams in food matrices, mostly for milk [11-14]. A number of methods have also been published for the analysis of β -lactams in tissue, including kidney [15-25], liver [20, 21, 24] and muscle [19-21, 23, 25-40]. The majority of these methods include a few penicillins and cephalosporins, or apply to specific drugs. The most comprehensive method was reported by Berendsen *et al.* for 22 β -lactams in poultry muscle, including penicillins, cephalosporins and carbapenems [2]. Becker *et al.* also developed a sensitive and specific method for 15 β -lactams in bovine milk, muscle and kidney [41]. However, in both cases, the sample preparation involves a laborious procedure with solid phase extraction (SPE) clean-up.

The objective of this study was to advance this analytical area and develop a comprehensive UHPLC-MS/MS method for monitoring the presence of the major β -lactams in both domestic and imported foods. The method aims to cover the antibiotics permitted for use in food-producing animals as listed under EU legislation [7], novel drugs that are considered as critically important antibiotics for human use only, and metabolites, such as desacetyl cephalirin (DAC), desfuroylceftiofur cysteine disulfide (DCCD) and desfuroylceftiofur dimer (DCD). A further objective of this study was to develop a fast and easy sample

preparation procedure that will reduce sample analysis time and simplify method implementation in regulatory laboratories.

2. Materials and Methods

2.1. Chemicals, standards and consumables

Ampicillin trihydrate, cloxacillin sodium salt monohydrate, dicloxacillin sodium salt hydrate, mecillinam, methicillin sodium salt, nafcillin sodium salt monohydrate, oxacillin sodium salt monohydrate, penicillin V potassium salt, cefalonium hydrate, cephapirin sodium, biapenem, doripenem monohydrate, meropenem trihydrate and faropenem sodium hydrate were purchased from Sigma-Aldrich (Dublin, Ireland). Amoxicillin trihydrate, penicillin G potassium salt, piperacillin, ticarcillin monosodium, cefadroxil hydrate, cefazolin sodium salt, cefoperazone sodium salt, cefotaxime sodium salt, cefquinome sulphate, cefuroxime sodium salt, cephalexin monohydrate and imipenem were purchased from LGC Standards (Teddington, Middlesex, UK). DCCD and DCD were a gift from Zoetis (Kalamazoo, MI, USA). Cefacetile, DAC sodium salt, ertapenem disodium 90% and the 10 internal standards amoxicillin-D₄, ampicillin-D₅, benzyl penicillanate-D₇ potassium salt (penicillin G-D₇), nafcillin-D₅ sodium salt, penicillin V-D₅, cefadroxil-D₄ (major), cefazolin-¹³C₂¹⁵N sodium salt, cephalexin-D₅ hydrate, DAC-D₆ sodium salt (major) and DCCD-D₃ were supplied by Toronto Research Chemicals (Toronto, ON, Canada).

Ultra-pure water (18.2 MΩ cm⁻¹) was generated in-house using a Millipore water purification system (Millipore, Cork, Ireland). Acetonitrile (MeCN), methanol (MeOH) and hexane (all UHPLC grade) were supplied by Romil Ltd (Cambridge, UK). Dimethyl sulfoxide (DMSO),

formic acid (HCOOH) 98-100%, trifluoroacetic acid (TFA), ammonium formate and ammonium acetate were supplied by Sigma-Aldrich. Glacial acetic acid (CH₃COOH) 100% was purchased from Merck (Darmstadt, Germany). Sorbents tested for d-SPE included NH₂ (Bondesil-NH₂ bulk sorbent, Agilent Technologies), PSA (Supelclean™ primary/secondary amine SPE bulk sorbent, Sigma-Aldrich), C₈ ISOLUTE bulk sorbent (Biotage, Uppsala, Sweeden), C₁₈ (Agilent Technologies Ltd, Cork, Ireland), Z-Sep, Z-Sep+ and Z-Sep/C₁₈ (Supel™ QuE, 15 mL tubes, Sigma-Aldrich). Anhydrous magnesium sulphate (MgSO₄) and sodium chloride (NaCl) were sourced from Sigma–Aldrich and Applichem (Darmstadt, Germany), respectively, and used as part of the QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) approach.

Polypropylene tubes (15 mL and 50 mL) were obtained from Sarstedt Ltd (Wexford, Ireland). Membrane filters (Captiva Econo Filter PTFE 13m 0.2 µm) were sourced from Agilent Technologies Ltd. A ME36S microbalance and an A200S digital electronic analytical balance, both from Sartorius (Dublin, Ireland), were used for standard preparation (weighing of standard material and required solvent, respectively). An Ultra-Turrax probe blender from IKA (Staufen, Germany), a MSE Mistral 3000i centrifuge from Davidson and Hardy (Dublin, Ireland) and a TurboVap LV evaporator from Biotage were used for sample preparation.

2.2. Preparation of standard solutions

Individual stock solutions of ampicillin, cloxacillin, dicloxacillin, mecillinam, methicillin, nafcillin, oxacillin, penicillin G, penicillin V, ticarcillin, cefadroxil, cephalixin, cephapirin, cefazolin, cefoperazone, cefotaxime, cefuroxime, biapenem, doripenem, ertapenem, imipenem, meropenem, faropenem and the internal standards ampicillin-D₅, penicillin G-D₇, penicillin V-D₅, nafcillin-D₅, cefadroxil-D₄, cefazolin-¹³C₂¹⁵N and cephalixin-D₅ were

prepared at a concentration of 1 mg mL⁻¹ in H₂O:MeCN (75:25, v/v). Stock solutions of piperacillin and cefalonium were prepared at a concentration of 1 mg mL⁻¹ in H₂O:MeCN (50:50, v/v). Amoxicillin, amoxicillin-D₄, DCCD and DCCD-D₃ were dissolved in H₂O:MeCN (50:50, v/v) at a concentration of 0.5 mg mL⁻¹. Cefacetile, cefquinome, DAC, DAC-D₆ and DCD were prepared in DMSO at a concentration of 1 mg mL⁻¹. After preparation, the solutions were stored in 2.5 mL aliquots in 15 mL polypropylene tubes at -80°C. Under these storage conditions, individual stock solutions were found to be stable for at least two months. Three intermediate standard solutions containing all the β-lactams at concentrations ranging from 4 µg mL⁻¹ to 200 µg mL⁻¹ were prepared in water. Working calibration solutions were prepared by diluting the intermediate stock solutions in water. A mixed internal standard solution was prepared at concentrations ranging from 0.2 µg mL⁻¹ to 10 µg mL⁻¹ in water. Working solutions were found to be stable for at least one month when stored at -80°C.

2.3. Sample preparation

Bovine muscle homogenate samples not containing any detectable β-lactam residues were used as negative controls. Sample aliquots (2 g ± 0.01 g) were weighed into a 50 mL polypropylene centrifuge tube. Extracted matrix calibrants were fortified with 100 µL of the working standard solutions. A 100 µL volume of internal standard solution was added to all calibrants, controls and test samples, which were then allowed to stand for 15 min. Water (1.9 mL for the extracted matrix calibrants and 2 mL for the controls and test samples) and MeCN (8 mL) were added in order to extract the analytes from the matrix, and the polypropylene tubes containing the samples were subsequently homogenised for 20 s over ice using an Ultra-Turrax probe blender. The homogenised samples were centrifuged for 15 min at 2842×g (4°C) and the supernatant was decanted into a 50 mL polypropylene tube containing

500 mg of C₁₈ sorbent. The samples were vortexed (40 s) and centrifuged for 15 min at 2842×g (4°C). A 5 mL aliquot of the supernatant was placed into a 15 mL polypropylene tube and evaporated under nitrogen on a TurboVap at 40°C to a final aqueous volume of < 1 mL. Solvent evaporation conditions required careful control to avoid degradation or loss of analytes. An evaporation temperature of 40°C was found to be optimal, while nitrogen gas pressure was dynamically controlled during the 55 min evaporation process. To avoid loss of sample, the nitrogen gas was maintained at 15 psi for the first 20 min. Subsequently, the pressure was gradually increased to 20 psi over 10 min, and the solvent evaporation time was greatly reduced by adapting this approach. After the evaporation, the volume was made up to 1 mL with water and the extracts were vortexed for 10 s prior to filtration through 0.2 µm PTFE syringe filters. The filtered extracts were collected directly into autosampler vials and a 10 µL volume was injected into the UHPLC-MS/MS system.

2.4. Instruments and UHPLC-MS/MS conditions

Samples were analysed using a Waters Acquity UHPLC system coupled to a Waters Quattro Premier triple quadrupole mass spectrometer (Milford MA, USA) equipped with ESI ionisation probe. The UHPLC-MS/MS system was controlled by MassLynx software (V4.1) and the results were processed by TargetLynx software (V4.1). Separation was performed on a stainless steel CSH C₁₈ analytical column (2.1 × 100 mm, particle size 1.7 µm) fitted with an in-line filter with a 0.2 µm pore size. The column was maintained at 30°C, as higher temperatures were found to negatively affect the precision of the method. A binary gradient separation comprising of 0.01% HCOOH and 0.2 mM ammonium acetate in H₂O (mobile phase A) and 0.01% HCOOH in MeCN (mobile phase B) was employed at a flow rate of 0.4 mL min⁻¹. The injection volume was 10 µL in full loop mode (overfill factor = 4). The gradient profile was as follows: (1) 0 - 1.5 min, 100% A; (2) 1.5 - 3.5 min, 80% A; (3) 3.5 -

8.5 min, 20% A; (4) 8.5 - 8.6 min, 0% A; (5) 8.6 - 10.5 min, 0% A; (6) 10.5 - 10.6 min, 100% A, which was held for 7.4 min (re-equilibration of the column), for a total run time of 18 min. The UHPLC autosampler was rinsed after each injection using strong (H₂O:MeCN, 20:80 v/v, 750 µL) and weak (0.01% HCOOH in H₂O:MeCN, 90:10 v/v, 1000 µL) washes. Sample temperature was maintained at 7°C in the autosampler. A divert valve was used to reduce source contamination (solvent delay: [a] 0 - 1.97 min; [b] 7.60 - 18 min). The mass spectrometer was operated in positive electrospray ionisation mode (ESI(+)) and the capillary voltage was set at 2.4 kV. Source and desolvation temperatures were 140°C and 450°C, respectively. Nitrogen was used as desolvation and cone gases and set at 1000 L h⁻¹ and 100 L h⁻¹, respectively. Argon collision gas flow rate was set at 0.22 mL min⁻¹. The MS tuning was performed by teed infusion of 1 µg mL⁻¹ standards/internal standards with mobile phase A:B (50:50, v/v). The cone voltage and collision energy were optimised for each analyte and two product ions were selected so that a minimum of four identification points were obtained for all the analytes as required by Commission Decision 2002/657/EC [42]. A multiple-reaction monitoring (MRM) method was developed with 15 different time-sectored events (Table 1). Inter-scan delay and inter-channel delay were set to 0.005 s.

2.5. Validation

The validation was performed following the 2002/657/EC guidelines [42]. The following parameters were assessed: identification, selectivity, linearity, trueness, within-laboratory repeatability (WLR), within-laboratory reproducibility (WLR), absolute recovery, decision limit (CC_α) and detection capability (CC_β). Matrix effects, limit of detection (LOD) and limit of quantitation (LOQ) were also evaluated as part of the validation process.

248 Identification was assessed by examining retention times, ion ratios and identification points.
249 The selectivity of the method was investigated through injecting standard solutions of all
250 analytes and internal standards individually and through testing 27 bovine muscles from
251 different animals, in order to check the presence of any interferences eluted at the retention
252 times of the analytes. The linearity of the curves was considered satisfactory if $R^2 \geq 0.98$ and
253 if individual residuals did not deviate by more than $\pm 20\%$ from the calibration curve. For the
254 MRL-substances, trueness, W_{Lr} and W_{LR} were assessed at 0.5, 1 and 1.5 times the MRLs
255 established by current legislation [7], except for DCD and DCCD, which were validated at
256 250, 500 and 750 $\mu\text{g kg}^{-1}$. The non-MRL substances were validated around a target level (TL)
257 that was identified based on the sensitivity of the method (Table 4). The W_{Lr} and W_{LR}
258 studies were both performed on three separate days by three different analysts. To evaluate
259 W_{Lr}, a different bovine muscle negative control was used on each day, and eight portions of
260 the sample were fortified at each validation level. To evaluate W_{LR}, 21 different bovine
261 muscle samples were used and seven bovine muscles, plus a mixture of them all, were
262 fortified at each validation level. The available labelled compounds were used in the
263 quantification of their corresponding analytes. Additionally, cloxacillin, dicloxacillin and
264 oxacillin were all corrected using nafcillin-D₅; cefacetile, cefotaxime, cefuroxime and
265 faropenem were corrected using cefazolin-¹³C₂¹⁵N; mecillinam and ertapenem using
266 ampicillin-D₅; cefalonium was corrected using cephalixin-D₅ and doripenem using DAC-D₆.
267 No internal standards were used for the other analytes, as the available labelled compounds
268 were found to be unsuitable in their quantification. Absolute recoveries were calculated by
269 comparing results from fortified samples to those of negative samples spiked post-extraction
270 at $2 \times$ lowest calibration level and $0.875 \times$ highest calibration level on three different
271 occasions. For the MRL-substances, CC α and CC β were calculated from the within-
272 laboratory reproducibility as defined in 2002/657/EC [42]. For the non-MRL substances, the

same approach was used considering the TL, in order to obtain clearly detectable $CC\alpha$ and $CC\beta$ values in the range of the calibration curves.

Matrix effects were also evaluated: 27 blank samples were spiked post-extraction at the MRLs or TLs and the signal obtained from those samples was compared to the signal obtained from a standard solution at the same concentration [43]. LODs and LOQs of all analytes were estimated from the blank samples fortified at the lowest calibration level on three different occasions, and measuring the signal-to-noise ratios at 3 and 10, respectively, for both quantifier and qualifier ions.

3. Results and Discussion

3.1. Method development

3.1.1. UHPLC-MS/MS conditions

The MS tuning experiments performed showed protonated molecules $[M+H]^+$ for most of the β -lactams. The monitored precursor ion for DCD was $[M+2H]^{2+}$; faropenem formed the sodium adduct $[M+Na]^+$, while cefacetrile and cefuroxime surprisingly formed the ammonium adduct $[M+NH_4]^+$. For the penicillins and cephalosporins, the product ions selected following low energy collision induced dissociation (CID) experiments generally matched those reported in the literature or characterized by Geis-Asteggianti *et al.* using electrospray ionization and quadrupole time-of-flight mass spectrometry (ESI-Q-TOF-MS) [44]. For the carbapenems, the monitored transitions were in agreement with those reported in the published LC-MS/MS methods for the determination of these drugs in urine, serum and

other biological samples [45-48]. In order to get optimal sensitivity, the following MS parameters were optimised: capillary voltage, cone voltage, collision energy, source temperature and desolvation temperature. The response for mecillinam and nafcillin was very intense, consequently the MS was detuned to avoid saturation of the signal at high concentrations.

The following additives were evaluated at different concentrations in both mobile phases A and B: 0.005%, 0.01% and 0.1% HCOOH; 0.01% and 0.1% CH₃COOH; 0.001% and 0.1% TFA. The addition of ammonium acetate and ammonium formate in mobile phase A only was also investigated at concentrations ranging from 0.1 to 0.5 mM. Also, three different UHPLC column chemistries were assessed, namely, BEH C₁₈, CSH C₁₈ and HSS T3. Optimal results were achieved when using a binary gradient separation with a mobile phase containing (A) 0.01% HCOOH and 0.2 mM ammonium acetate in H₂O and (B) 0.01% HCOOH in MeCN, on a CSH C₁₈ column. A concentration of 0.01% HCOOH provided the best overall response, while higher concentrations of the same acid or use of different acids decreased sensitivity. Moreover, the removal of TFA from the UHPLC system after its use was found to be very difficult, and its presence seemed to affect other assays performed on the same instrument. The response for all analytes was also decreased when using higher concentrations of salts, consequently 0.2 mM ammonium acetate was added to mobile phase A to promote the formation of the ammonium adduct for two cephalosporins, namely cefacetile and cefuroxime. The results obtained were in agreement with what reported in the literature. Formic acid is, indeed, the most common acidic agent used in the analysis of β -lactams [11-13]. In addition, Becker *et al.* evaluated the effect of different concentrations of HCOOH and observed that the optimal chromatographic separation and ionisation efficiency were obtained with small amounts of acid in the eluent [41].

323

324 Different column temperatures were also tested (30°C, 40°C and 50°C) and the lower
325 temperature was chosen, as it was shown to improve the precision of the method. The
326 analytes particularly affected by the column temperature were imipenem, faropenem,
327 cefuroxime, DCCD and DCD. A series of 25 injections of a 50 ng mL⁻¹ mixed standard
328 solution showed that the relative standard deviation (RSD%) for these compounds ranged
329 from 6.3% (cefuroxime) to 11.5% (imipenem) when the column temperature was set at 50°C,
330 while the precision improved at 30°C, with RSDs% ranging from 1.4% (faropenem) to 4.0%
331 (DCD). The autosampler temperature was set at 7°C to prevent a potential degradation of the
332 analytes. During each run performed in the proposed study, a matrix sample fortified at the
333 MRLs/TLs was injected as check sample at evenly spaced intervals throughout the run to
334 monitor the stability of the analytes in the autosampler. All the compounds were found to be
335 stable for at least 24 hours at 7°C, as the absolute area obtained from the check sample stayed
336 constant until the end of the run and the RSD% was not greater than 7.9%, except for
337 imipenem (RSD% = 10.2%), for which the peak area decreased gradually, due to the
338 instability of this compound.

339

340 Gradient conditions and flow rate were adjusted in order to reach optimal chromatographic
341 separation. The multiple reaction monitoring chromatograms in Fig. 2 show that all the 30 β-
342 lactams and the 10 internal standards were separated within the first 8 min of the
343 chromatographic run. However, after 8.6 min, the gradient profile was held at 100% B for
344 approximately 2 min, in order to remove the most non-polar matrix components from the
345 column. Subsequently, a 7.4-min hold at 100% A was necessary for column equilibration,
346 resulting in a total cycle time of 18 min.

347

3.1.2. Sample preparation

One of the main objectives of this study was to develop a rapid and easy sample preparation procedure. The QuEChERS method is well known for being a streamlined approach and was evaluated during preliminary studies, but gave unacceptably low recovery of polar analytes, likely due to partitioning into the aqueous layer. In contrast, a suitable protocol was identified based on the methods developed by Fagerquist *et al.* [17] and Mastovska and Lightfield [22]. They found that an extraction solvent composed of MeCN:H₂O (80:20, v/v) was optimal for the isolation of 11 β -lactams from 1 g of bovine kidney. The sample weight was subsequently increased to 2 g by Schneider *et al.* [49] and the extraction procedure was applied to bovine kidney for the analysis of 120 veterinary drugs from different classes, including seven penicillins and eight cephalosporins. Geis-Asteggianti *et al.* [50] also adapted the same sample preparation to 2 g of bovine muscle for monitoring >100 veterinary drug residues, with the inclusion of seven penicillins and three cephalosporins. In our method, the volume of the extraction solvent was kept at 10 mL and the procedure used for the selective and efficient isolation of 12 penicillins, 12 cephalosporins, five carbapenems and faropenem from 2 g of bovine muscle samples.

With the exception of cefazolin, cefotaxime and faropenem, for which no significant difference was observed, all the analytes seemed to be affected by the sample temperature, which needed to be carefully controlled along the sample preparation process to prevent degradation of the β -lactam residues. During the method development, it was found that, after the addition of the extraction solvent, samples required homogenisation over ice, and that the maximum solvent evaporation temperature had to be restricted to 40°C.

372 A number of sorbent materials were investigated for d-SPE: NH₂ (500 mg), PSA (500 mg),
373 C₁₈ (500 mg), C₈ (500 mg) and zirconium dioxide-based sorbents (Z-Sep (500 mg); Z-Sep+
374 (500 mg); Z-Sep/C₁₈ (120 mg/ 300 mg)). NH₂ and PSA are sorbents effective at retaining
375 fatty acids and other organic acids present in food. C₁₈ is often used for removing lipophilic
376 compounds in food of animal origin, which has high lipid content. C₈ is a sorbent generally
377 employed for the extraction of non-polar compounds from matrices with low-fat content. Z-
378 Sep is a novel bonded silica sorbent that has been developed for the selective isolation of
379 hydrophobic analytes from fatty matrices, while Z-sep/C₁₈ or Z-sep+ are recommended for
380 samples containing <15% fat and >15% fat, respectively. A liquid-liquid partitioning clean-
381 up was also tested using 6 mL of hexane, as an alternative to the d-SPE clean-up step. Geis-
382 Asteggianti *et al.* [50] reported a similar study for several veterinary drugs in bovine muscle,
383 including 11 β -lactams. In their work, they evaluated the efficiency of different d-SPE
384 sorbents (C₁₈, Z-Sep, Z-Sep+) and/or partitioning with hexane in removing matrix co-
385 extracting by gravimetric measurements. The liquid-liquid partitioning with hexane was
386 found to be the least efficient, followed by C₁₈ clean-up. A combination of Z-Sep and C₁₈
387 resulted in the highest matrix removal, followed by Z-Sep + C₁₈ + hexane and Z-Sep alone.
388 However, further investigation on matrix effects and recoveries discouraged the use of these
389 clean-up strategies; consequently, a combination of C₁₈ and hexane was used in their sample
390 preparation. In the proposed study, sensitivity, absolute recoveries and precision were
391 evaluated for each of the clean-up approaches tested. NH₂ and PSA sorbents were
392 immediately discarded, as they provided unacceptable results for the majority of the analytes.
393 Mastovska and Lightfield (2008) also found PSA to be unsuitable for the analysis of β -
394 lactams in bovine kidney, finding these drugs were retained by PSA, likely due to the
395 carboxylic group present in the molecular structure [22]. Hexane provided high recoveries for
396 cloxacillin, nafcillin and oxacillin, while the recoveries were not satisfactory for cefuroxime

and DAC when using C₈ sorbent. Among the zirconium dioxide-based sorbents, Z-Sep showed the lowest overall recoveries (particularly for the carbapenems) (Fig. 3). This result was in good agreement with that reported by Geis-Asteggianti *et al.* [50]. In their study, they found that the β -lactams were strongly retained when using a combination of Z-Sep and hexane. This could be due to the fact that zirconium can easily react with the carboxyl and hydroxyl groups present in the β -lactam structures and form strong bonds. In the proposed work, C₁₈, Z-Sep+ and Z-Sep/C₁₈ were found to be the most efficient strategies for all the substances, except for imipenem, which was found to be retained also by Z-Sep+ and Z-Sep/C₁₈ (Fig. 4). In addition, more precise results were obtained for the majority of compounds when using the C₁₈ sorbent. As a consequence, C₁₈ was chosen for the proposed sample preparation.

3.2. Standard stability

The stability issues associated with penicillin methanolic stock solutions were discovered when tuning was performed. The standards were initially dissolved in MeOH or H₂O/MeOH and stored at -30°C. The MS tuning of 1 $\mu\text{g mL}^{-1}$ standard aqueous solutions was originally performed immediately after standard preparation and subsequently repeated two weeks later. Interestingly, the MS spectra obtained from the second experiment for cloxacillin, dicloxacillin and penicillin V showed a mass corresponding to their methyl ester $[\text{M}+\text{CH}_3\text{OH}+\text{H}]^+$. This was likely due to a breakdown of the β -lactam ring and the formation of a methanolic degradation product. The results raised curiosity, mostly because MeOH is a common solvent used in the analysis of β -lactams. In addition, contradictory information had been found in the literature regarding the stability of penicillin standards in solvent [17, 22, 41, 51-53]. It was decided, then, to investigate the phenomenon for one of the problematic analytes, i.e. cloxacillin. The MS tuning was repeated after preparation of 1 mg mL^{-1} stock

solutions in MeOH and H₂O:MeCN (75:25, v/v), which were stored for 11 days in 1 mL aliquots under different storage conditions (room temperature, +4°C, -30°C and -80°C). The spectra obtained from the stock solutions in MeOH showed a mass at 436.0 m/z ([M+H]⁺), but also a mass at 468.0 m/z ([M+CH₃OH+H]⁺) that did not appear in the spectra obtained from the standard dissolved in H₂O:MeCN (Fig. 5). All the different aliquots were subsequently diluted in water at a concentration of 10 µg mL⁻¹. A MS scan was acquired following column injections and the chromatograms extracted for the ions 436.0 and 468.0 m/z . As shown in Fig. 6, the signal obtained from 436.0 m/z was of lower intensity for the standard prepared in MeOH compared to the one prepared in H₂O:MeCN, indicating some degradation for cloxacillin. The degradation in aliquots that had been stored at lower temperature was less (22% degradation at -80°C and 19% at -30°C) than the aliquots stored at +4°C (81% degradation) or room temperature (99% degradation). However, a fresh solution prepared in MeOH also showed a mass at 468.0 m/z and 31% degradation for 436.0 m/z , which was surprisingly higher than what observed after storage at -80°C and -30°C. This was likely due to the fact that the fresh solution had been left at room temperature for approximately three hours before the dilution in water and the column injection. It can be concluded that the degradation of cloxacillin starts as soon as the standard is dissolved in MeOH and seems to be faster at higher storage temperatures. It should be noted that no significant difference was observed between the aliquots stored at +4°C and room temperature, but this was likely due to the saturation of the detector after the injection of high concentration standard solutions.

The results obtained from the experiment were in good agreement with what was reported by Tyczkowska *et al.* [53] and Pellicciotti *et al.* [54]. Although no degradation appeared to be observed for the carbapenems and the cephalosporins, MeOH was completely removed from

the standard preparation procedure, in order to avoid the presence of this solvent in the mixed standard solutions. As a consequence, alternative solvents (H₂O:MeCN or DMSO) were selected. In addition, the standard solutions were stored in small aliquots (2.5 mL) at -80°C.

3.3. Method validation

3.3.1. Identification and confirmatory criteria

In order to satisfy the confirmatory criteria outlined in 2002/657/EC, three identification points (one point for the precursor ion and 1.5 points for each product ion) are required to identify a group B substance. This criteria was exceeded in the proposed method, with four identification points obtained for each analyte. All the retention times and ion ratios measured during the validation study were within the maximum permitted tolerances [42], showing the suitability of the method in the identification and confirmation of 30 β -lactams in bovine muscle samples.

3.3.2. Selectivity, linearity, LODs and LOQs

Analytes and internal standards were injected individually in order to monitor for interferences in the UHPLC-MS/MS traces. The effects of cross-talk and isobaric interferences were minimised by identifying unique product ions for each compound and optimising the chromatographic separation, respectively. During the selectivity study, no matrix interferent peaks were observed for the majority of the analytes. An interference was observed for cefazolin, cefadroxil and cefadroxil-D₄ product ions at 323.1, 114.1 and 213.1 m/z , respectively. As a consequence, alternative product ions of lower intensity were selected. An interference was also observed for doripenem and cefalonium at 274.2 and 337.0 m/z , respectively. In the case of cefalonium, only two product ions (152.0 and 337.0 m/z) were

generated in low energy CID experiments. In contrast, doripenem gave two additional product ions (112.1 and 318.1 m/z) to those selected (342.2 and 274.2 m/z), but the choice of these ions would have significantly decreased the sensitivity of the method. In addition, the quantifier and the qualifier ions for both analytes were of similar intensity; consequently it could be concluded that the observed peak was a real matrix interference (and not the actual analyte coming from a drug treatment of the animals), as no peaks were observed for the second ion. It is highlighted that the absolute area count obtained for the interferences was found to be negligible in all the tested blank samples, however further investigation would be needed for doripenem and cefalonium by analysing a larger number of bovine muscle samples.

Linearity was achieved over the calibration range of the method. R^2 was ≥ 0.988 for all the analytes, while individual residuals were in the $\pm 20\%$ range from the calibration curve. LODs ranged between 0.010 $\mu\text{g kg}^{-1}$ and 4.8 $\mu\text{g kg}^{-1}$; LOQs ranged from 0.035 $\mu\text{g kg}^{-1}$ to 16.0 $\mu\text{g kg}^{-1}$ (Table 2).

3.3.3. Matrix effects

During the matrix effect study, ion suppression was observed for the majority of the analytes, while biapenem, cefoperazone, cefquinome, imipenem and methicillin showed ion enhancement. The greatest amount of ion suppression and enhancement was observed for amoxicillin and methicillin, respectively (75% and 73%). The use of the 10 internal standards significantly improved the precision of the method (Table 2 and Table 3) and represented a distinguishing feature of the proposed work. With the exception of Berendsen *et al.* [2], none of the previously published methods included such a large number of deuterated standards. However, the available labelled compounds were found not suitable in the quantification of

biapenem, cefoperazone, cefquinome, DCD, imipenem, meropenem, methicillin, piperacillin and ticarcillin. Among the analytes corrected using a labelled compound, mecillinam showed the greatest variability (RSD% = 26.2%). It is highlighted that no MRL is required for mecillinam in bovine species. Cefquinome was the only MRL-substance that showed ion enhancement. The use of an internal standard is particularly important for those analytes for which ion enhancement is observed. In the event of positive samples being found, this matrix effect can potentially lead, indeed, to false non-compliant results. However, in the proposed method, the available labelled compounds were found to be unsuitable for the quantification of cefquinome. In addition, the matrix effect study showed a variability of 28.9%. This problem should be addressed in the future through the implementation of a suitable isotopically labelled standard.

3.3.4. Trueness, precision, absolute recovery, CCa and CCβ

The 2002/657/EC guidelines specify that the trueness for mass fractions between 1 and 10 µg kg⁻¹ should be in the range 70-110%, while the trueness for mass fractions ≥ 10 µg kg⁻¹ should be in the range 80-110%. The trueness was satisfactory for all the analytes under WLR conditions (range of 86-108%). On day two of the WLR study, the mixed internal standard solution was accidentally added twice to one of the samples fortified at 1.5 MRL/TL, while the internal standard was not added to the following sample. Consequently, six replicates were used on that day for the assessment of WLR at 1.5 MRL/TL for those analytes corrected using an internal standard. The trueness was satisfactory for the majority of substances under WLR conditions (range of 81-108%). The exceptions were mecillinam (69% accuracy at TL), meropenem (115-125%) and faropenem (119-143%). For mecillinam and faropenem, the trueness was outside the acceptable range even though the internal standards cefazolin-¹³C₂¹⁵N and ampicillin-D₅ were used, respectively. In the case of meropenem, for which no

internal standards were used, the corresponding labelled compound, namely meropenem-D₆, could be implemented into the method to address the problem. The 2002/657/EC guidelines also state that the coefficient of variations (CVs) for mass fractions between 100 and 1000 µg kg⁻¹ shall be ≤ 23% in WLR conditions. For mass fractions lower than 100 µg kg⁻¹, the CVs shall be as low as possible. In addition, the CVs calculated under WLR conditions should not be greater than two thirds of the CVs calculated under WLR conditions. The precision of the method was satisfactory for all analytes in both WLR (2.0-29.9%) and WLR studies (1.5-17.3%), as the CVs did not exceed the maximum CVs calculated by the Horwitz equation and two thirds of these values, respectively.

CCα ranged from 10.6 µg kg⁻¹ (cefazolin) to 677 µg kg⁻¹ (DCD), depending on the analyte. Overall absolute recoveries (OAR) ranged from 61% to 89% for all the compounds, except for cefadroxil, DCCD, DCD, imipenem and ertapenem (Table 3). Cefadroxil, ertapenem and imipenem are relatively polar analytes, and the lower recovery rates were probably due to the incomplete extraction from the matrix, caused by the specific extraction solvent composition used in the sample preparation. In the case of DCCD and DCD, the disulphide bonds present in their molecular structures likely undergo exchange with protein thiols or other disulphide bonds in the homogenised tissue samples, resulting in losses during the deproteinization step or in the non-detection by the LC-MS/MS system [17]. It is important to state that all active metabolites should be included in the analysis of ceftiofur residues in food products, as the MRL for muscle tissue in all mammalian food-producing species is defined as the sum of all residues retaining the β-lactam structure, expressed as desfuroylceftiofur [7]. After intramuscular administration, ceftiofur is rapidly metabolised to desfuroylceftiofur and further products, including protein-bound desfuroylceftiofur [26]. Since no hydrolysis is involved, the method described in this study does not detect protein-bound metabolites, consequently it leads to an underestimation of the total ceftiofur content. To measure the total

amount of ceftiofur present in a sample, all the disulphide bonds should be reduced and desfuroylceftiofur released, followed by stabilisation of the thiol group by derivatization [22, 26].

Cephapirin could not be directly included in the proposed method. A very rapid conversion of cephapirin to its active metabolite (DAC) in milk, kidney homogenate and fortified muscle tissue samples was reported in literature [16, 41, 55, 56]. The preliminary experiments performed during the development of the proposed method employed the use of mixes containing both cephapirin and DAC and showed no absolute recovery for cephapirin and 180-190% absolute recovery for DAC. For this reason, a quantitative conversion of cephapirin to DAC was hypothesised after fortification of blank muscle tissue homogenate samples (in approximately 15 min), resulting in the impossibility to quantify cephapirin directly. However, as highlighted by Fagerquist and Lightfield (2003), if this compound was present in a real sample, it would be rapidly and completely converted to DAC after a preliminary homogenisation aimed to obtain a representative sample [16].

Although no MRLs are required for cefacetrile, cefalonium, cefazolin, cefoperazone, mecillinam and penicillin V in bovine muscle, the inclusion of these drugs in the proposed method could provide useful information on their transmission through the food chain in the case that positive samples are found. In addition, as a MRL is required for penicillin V in liver, kidney, skin and fat, and for cefacetrile, cefalonium, cefazolin and cefoperazone in milk, the sample preparation could be adapted to different matrices, while using the same UHPLC-MS/MS method.

4. Comparison with other existing methods

Relatively few multi-residue LC-MS methods have been published for the analysis of β -lactams in muscle tissue. A preliminary literature review found few methods dealing with penicillins and cephalosporins [19-21, 25, 27-40]. Most of these methods include few analytes or exclude cephapirin and ceftiofur metabolites. On the other hand, other methods describe dedicated tests for ceftiofur and its metabolites only [23, 26]. With the exception of Berendsen *et al.* [2], no methods have been reported at all for carbapenems in food products. Although the Berendsen method represents the most valid approach for monitoring ceftiofur and the total amount of its metabolites, it involves a complex derivatisation process followed by SPE clean-up. A time consuming sample preparation including SPE clean-up is also involved in the method reported by Becker *et al.* [41] (Fig. 7). The method described in this study consists of a fast and easy sample preparation protocol and allows the simultaneous analysis of 30 β -lactam antibiotics, including carbapenems, DAC and two of the ceftiofur metabolites, in a reasonably short chromatographic run (18 min). The method published by Geis-Asteggianti *et al.* (2012) had previously applied the same sample preparation to bovine muscle for the analysis of several veterinary drug residues belonging to different classes of antibiotics [50]. It has to be said that the development of a multi-residue method is particularly difficult when several analytes with a wide range of chemical properties are present. However, for the 10 β -lactams included in the Geis-Asteggianti method, the average reproducibility %CV ranged from 8.0% to 49.0%, while the reproducibility of the proposed method could be improved by using 10 different internal standards and ranged between 2.0% and 29.9%. Finally, limits of quantitation for the MRL-substances were well below half the MRLs set by current legislation [7].

5. Conclusions

This research has demonstrated that it is possible to accurately measure 30 β -lactam residues belonging to the penicillin, cephalosporin and carbapenem families in bovine muscle samples by UHPLC-MS/MS following a generic sample preparation protocol. A range of different clean-up procedures were evaluated during method development, and d-SPE with C₁₈ was found to give the best overall recoveries. The performance of the method was assessed during matrix effects studies and validation studies, showing it was necessary to include isotopically labelled internal standards to improve accuracy. It was shown in this work that β -lactam residues are particularly unstable and protocols need to be implemented in test methods to mitigate against this problem. Firstly, this work demonstrated that MeOH should not be employed for the preparation of standards, as it causes degradation of some penicillins. This can potentially lead to the identification and selection of incorrect precursor/product ions, as well as an overestimation of beta-lactam residue levels in test samples. This standard stability problem can be addressed through the use of alternative solvents, such as H₂O/MeCN and DMSO. Additionally, it is recommended that MeOH should not be used during sample preparation or during UHPLC-MS/MS analysis. This research concurs with previous research findings that proposed that temperature needs to be carefully controlled during β -lactam analysis. Our study found that analyte recovery and method precision can be improved by extraction over ice, maintaining evaporation at 40°C and reducing chromatographic column temperature to 30°C.

The final method is advantageous over previous methods because it includes the widest range of β -lactams and allows the analysis of carbapenem residues without the need for derivatisation. The method has a reasonably high throughput that allows one analyst working alone to analyse 30 test samples along with calibration and quality control samples in one

day. The method was validated according to the 2002/657/EC guidelines and was demonstrated to be fit for purpose.

Acknowledgements

This research was funded by the FIRM programme administered by the Irish Department of Agriculture, Food and the Marine (Contract 13/F484). We would like to thank Simona Pellicciotti and Guglielmo Dusi for their advice on the stability of standard solutions (Istituto Zooprofilattico della Lombardia e dell'Emilia Romagna, Brescia, Italy).

References

- [1] P.T. Reeves, Antibiotics: Groups and Properties, in: J. Wang, J.D. MacNeil, J.F. Kay (Eds.), Chemical Analysis of Antibiotic Residues in Food, Wiley (2012) pp. 1-60.
- [2] B.J. Berendsen, H.W. Gerritsen, R.S. Wegh, S. Lameris, R. van Seville, A.A. Stolker, M.W. Nielen, Comprehensive analysis of β -lactam antibiotics including penicillins, cephalosporins, and carbapenems in poultry muscle using liquid chromatography coupled to tandem mass spectrometry, Analytical and bioanalytical chemistry 405(24) (2013) 7859-74.
- [3] F.C. Tenover, Mechanisms of Antimicrobial Resistance in Bacteria, The American Journal of Medicine 119(6, Supplement 1) (2006) S3-S10.
- [4] R.E. Hornish, S.F. Kotarski, Cephalosporins in veterinary medicine - ceftiofur use in food animals, Curr Top Med Chem 2(7) (2002) 717-31.
- [5] R.N. Jones, Resistance Patterns Among Nosocomial Pathogens: Trends Over the Past Few Years, Chest 119(2, Supplement) (2001) 397S-404S.

645 [6] Y. Pfeifer, A. Cullik, W. Witte, Resistance to cephalosporins and carbapenems in Gram-
646 negative bacterial pathogens, *International Journal of Medical Microbiology* 300(6) (2010)
647 371-379.

648 [7] Commission Regulation (EU) No 37/2010 of 22 December 2009 on pharmacologically
649 active substances and their classification regarding maximum residue limits in food-stuffs of
650 animal origin. *Official Journal of the European Union* Nr. 15, 20.1.2010 (2010).

651 [8] Commission Notice (2015/C 299/04) on guidelines for the prudent use of antimicrobials
652 in veterinary medicine. *Official Journal of the European Union* Nr. 299, 11.9.2015 (2015).

653 [9] K.M. Papp-Wallace, A. Endimiani, M.A. Taracila, R.A. Bonomo, Carbapenems: past,
654 present, and future, *Antimicrobial agents and chemotherapy* 55(11) (2011) 4943-60.

655 [10] P. Gowik, S. Uhlig, K. Frost, A monitoring study on the level and frequency of antibiotic
656 residues in Germany, *Conference on Residues of Veterinary Drugs in Food (EuroResidue*
657 *VIII)*, Egmond aan Zee (The Netherlands), 23-25 May 2016.

658 [11] A. Gentili, D. Perret, S. Marchese, Liquid chromatography-tandem mass spectrometry
659 for performing confirmatory analysis of veterinary drugs in animal-food products, *Trends in*
660 *Analytical Chemistry* 24(7) (2005) 704-733.

661 [12] L. Kantiani, M. Farre, D. Barcelo, Analytical methodologies for the detection of beta-
662 lactam antibiotics in milk and feed samples, *Trends in Analytical Chemistry* 28(6) (2009)
663 729-744.

664 [13] F.J. Lara, M. del Olmo-Iruela, C. Cruces-Blanco, C. Quesada-Molina, A.M. Garcia-
665 Campana, Advances in the determination of beta-lactam antibiotics by liquid
666 chromatography, *Trends in Analytical Chemistry* 38 (2012) 52-66.

667 [14] V. Samanidou, S. Nisyrliou, Multi-residue methods for confirmatory determination of
668 antibiotics in milk, *Journal of separation science* 31(11) (2008) 2068-2090.

669 [15] M.P. de Almeida, C.P. Rezende, F.D. Ferreira, L.F. de Souza, D.Cristina S. de Assis,
670 T.C. de Figueiredo, M. de Oliveira Leite, S. de Vasconcelos Cançado, Optimization and
671 validation method to evaluate the residues of β -lactams and tetracyclines in kidney tissue by
672 UPLC–MS/MS, *Talanta* 144 (2015) 922-932.

673 [16] C.K. Fagerquist, A.R. Lightfield, Confirmatory analysis of beta-lactam antibiotics in
674 kidney tissue by liquid chromatography/electrospray ionization selective reaction monitoring
675 ion trap tandem mass spectrometry, *Rapid Communications in Mass Spectrometry* 17(7)
676 (2003) 660-671.

677 [17] C.K. Fagerquist, A.R. Lightfield, S.J. Lehotay, Confirmatory and quantitative analysis of
678 beta-lactam antibiotics in bovine kidney tissue by dispersive solid-phase extraction and liquid
679 chromatography-tandem mass spectrometry, *Anal. Chem.* 77(5) (2005) 1473-1482.

680 [18] S. Feng, C. Chattopadhyaya, P. Kijak, O.A. Chiesa, E.A. Tall, A determinative and
681 confirmatory method for ceftiofur metabolite desfuroylceftiofur cysteine disulfide in bovine
682 kidney by LC-MS/MS, *J Chromatogr B* 898 (2012) 62-8.

683 [19] K. Granelli, C. Branzell, Rapid multi-residue screening of antibiotics in muscle and
684 kidney by liquid chromatography-electrospray ionization-tandem mass spectrometry,
685 *Analytica Chimica Acta* 586(1-2) (2007) 289-95.

686 [20] Y. Ito, Y. Ikai, H. Oka, H. Matsumoto, Y. Miyazaki, K. Takeba, H. Nagase, Application
687 of ion-exchange cartridge clean-up in food analysis IV. Confirmatory assay of
688 benzylpenicillin, phenoxymethylpenicillin, oxacillin, cloxacillin, nafcillin and dicloxacillin in
689 bovine tissues by liquid chromatography-electrospray ionization tandem mass spectrometry,
690 *Journal of chromatography. A* 911(2) (2001) 217-223.

691 [21] Y. Ito, T. Goto, H. Oka, H. Matsumoto, K. Takeba, Application of ion-exchange
692 cartridge clean-up in food analysis - VI. Determination of six penicillins in bovine tissues by

693 liquid chromatography-electrospray ionization tandem mass spectrometry, *Journal of*
694 *chromatography. A* 1042(1-2) (2004) 107-111.

695 [22] K. Mastovska, A.R. Lightfield, Streamlining methodology for the multiresidue analysis
696 of beta-lactam antibiotics in bovine kidney using liquid chromatography-tandem mass
697 spectrometry, *Journal of chromatography. A* 1202(2) (2008) 118-23.

698 [23] S. Mompelat, M.-P. Fourmond, M. Laurentie, E. Verdon, D. Hurtaud-Pessel, J.-P.
699 Abjean, Validation of a liquid chromatography–high-resolution mass spectrometry method
700 for the analysis of ceftiofur in poultry muscle, kidneys and plasma: A unique accuracy profile
701 for each and every matrix, *Journal of Chromatography A* 1407 (2015) 119-129.

702 [24] T.A.M. Msagati, M.M. Nindi, Determination of β -lactam residues in foodstuffs of
703 animal origin using supported liquid membrane extraction and liquid chromatography–mass
704 spectrometry, *Food chemistry* 100(2) (2007) 836-844.

705 [25] F. van Holthoon, P.P. Mulder, E.O. van Bennekom, H. Heskamp, T. Zuidema, H.J. van
706 Rhijn, Quantitative analysis of penicillins in porcine tissues, milk and animal feed using
707 derivatisation with piperidine and stable isotope dilution liquid chromatography tandem mass
708 spectrometry, *Analytical and bioanalytical chemistry* 396(8) (2010) 3027-40.

709 [26] B.J. Berendsen, L.A. Stolker, M.W. Nielen, Assessment of liquid chromatography-
710 tandem mass spectrometry approaches for the analysis of ceftiofur metabolites in poultry
711 muscle, *Food Additives & Contaminants: Part A* 29(2) (2012) 197-207.

712 [27] D. Chen, J. Yu, Y. Tao, Y. Pan, S. Xie, L. Huang, D. Peng, X. Wang, Y. Wang, Z. Liu,
713 Z. Yuan, Qualitative screening of veterinary anti-microbial agents in tissues, milk, and eggs
714 of food-producing animals using liquid chromatography coupled with tandem mass
715 spectrometry, *Journal of Chromatography B* 1017–1018 (2016) 82-88.

716 [28] C. ChiaoChan, U. Koesukwiwat, S. Yudthavorasit, N. Leepipatpiboon, Efficient
717 hydrophilic interaction liquid chromatography-tandem mass spectrometry for the multiclass

718 analysis of veterinary drugs in chicken muscle, *Analytica Chimica Acta* 682(1-2) (2010) 117-
719 29.

720 [29] J. Chico, A. Rubies, F. Centrich, R. Companyo, M.D. Prat, M. Granados, High-
721 throughput multiclass method for antibiotic residue analysis by liquid chromatography-
722 tandem mass spectrometry, *Journal of chromatography. A* 1213(2) (2008) 189-199.

723 [30] K. Granelli, C. Elgerud, A. Lundstrom, A. Ohlsson, P. Sjoberg, Rapid multi-residue
724 analysis of antibiotics in muscle by liquid chromatography-tandem mass spectrometry,
725 *Analytica Chimica Acta* 637(1-2) (2009) 87-91.

726 [31] M.P. Hermo, P. Gómez-Rodríguez, J. Barbosa, D. Barrón, Metabolomic assays of
727 amoxicillin, cephalixin and ceftiofur in chicken muscle: Application to treated chicken
728 samples by liquid chromatography quadrupole time-of-flight mass spectrometry, *Journal of*
729 *pharmaceutical and biomedical analysis* 85 (2013) 169-178.

730 [32] Z. Huang, X.-D. Pan, B.-f. Huang, J.-J. Xu, M.-L. Wang, Y.-P. Ren, Determination of 15
731 β -lactam antibiotics in pork muscle by matrix solid-phase dispersion extraction (MSPD) and
732 ultra-high pressure liquid chromatography tandem mass spectrometry, *Food Control* 66
733 (2016) 145-150.

734 [33] D. Hurtaud-Pessel, T. Jagadeshwar-Reddy, E. Verdon, Development of a new screening
735 method for the detection of antibiotic residues in muscle tissues using liquid chromatography
736 and high resolution mass spectrometry with a LC-LTQ-Orbitrap instrument, *Food Additives*
737 *& Contaminants: Part A* 28(10) (2011) 1340-51.

738 [34] W. Li, H. Shen, Y. Hong, Y. Zhang, F. Yuan, F. Zhang, Simultaneous determination of
739 22 cephalosporins drug residues in pork muscle using liquid chromatography–tandem mass
740 spectrometry, *Journal of Chromatography B* 1022 (2016) 298-307.

741 [35] C.A. Macarov, L. Tong, M. Martinez-Huelamo, M.P. Hermo, E. Chirila, Y.X. Wang, D.
742 Barron, J. Barbosa, Multi residue determination of the penicillins regulated by the European

743 Union, in bovine, porcine and chicken muscle, by LC-MS/MS, Food chemistry 135(4) (2012)
744 2612-21.

745 [36] S. Moretti, G. Dusi, D. Giusepponi, S. Pellicciotti, R. Rossi, G. Saluti, G. Cruciani, R.
746 Galarini, Screening and confirmatory method for multiclass determination of 62 antibiotics in
747 meat, Journal of Chromatography A 1429 (2016) 175-188.

748 [37] R. Perez-Burgos, E.M. Grzelak, G. Gokce, J. Saurina, J. Barbosa, D. Barron, Quechers
749 methodologies as an alternative to solid phase extraction (SPE) for the determination and
750 characterization of residues of cephalosporins in beef muscle using LC-MS/MS, J
751 Chromatogr B 899 (2012) 57-65.

752 [38] C.P. Rezende, M.P. Almeida, R.B. Brito, C.K. Nonaka, M.O. Leite, Optimisation and
753 validation of a quantitative and confirmatory LC-MS method for multi-residue analyses of
754 beta-lactam and tetracycline antibiotics in bovine muscle, Food Additives & Contaminants:
755 Part A 29(4) (2012) 541-9.

756 [39] Y.Y. Tang, H.F. Lu, H.Y. Lin, Y.C. Shin, D.F. Hwang, Development of a Quantitative
757 Multi-Class Method for 18 Antibiotics in Chicken, Pig, and Fish Muscle using UPLC-
758 MS/MS, Food Anal. Methods 5(6) (2012) 1459-1468.

759 [40] M.E. Dasenaki, C.S. Michali, N.S. Thomaidis, Analysis of 76 veterinary
760 pharmaceuticals from 13 classes including aminoglycosides in bovine muscle by hydrophilic
761 interaction liquid chromatography–tandem mass spectrometry, Journal of Chromatography A
762 1452 (2016) 67-80.

763 [41] M. Becker, E. Zittlau, M. Petz, Residue analysis of 15 penicillins and cephalosporins in
764 bovine muscle, kidney and milk by liquid chromatography–tandem mass spectrometry,
765 Analytica Chimica Acta 520(1–2) (2004) 19-32.

766 [42] Commission Decision 2002/657/EC of 12 August 2002 implementing Council Directive
767 96/23/EC concerning performance of analytical methods and the interpretation of results.
768 Official Journal of the European Union Nr. 221, 17.8.2002 (2002).

769 [43] H. Trufelli, P. Palma, G. Famiglini, A. Cappiello, An overview of matrix effects in
770 liquid chromatography–mass spectrometry, *Mass Spectrometry Reviews* 30(3) (2011) 491-
771 509.

772 [44] L. Geis-Asteggianti, A. Nunez, S.J. Lehotay, A.R. Lightfield, Structural characterization
773 of product ions by electrospray ionization and quadrupole time-of-flight mass spectrometry to
774 support regulatory analysis of veterinary drug residues in foods, *Rapid Commun Mass*
775 *Spectrom* 28(10) (2014) 1061-1081.

776 [45] R. Cazorla-Reyes, R. Romero-Gonzalez, A.G. Frenich, M.A. Rodriguez Maresca, J.L.
777 Martinez Vidal, Simultaneous analysis of antibiotics in biological samples by ultra high
778 performance liquid chromatography-tandem mass spectrometry, *Journal of pharmaceutical*
779 *and biomedical analysis* 89 (2014) 203-12.

780 [46] S. Lefevre, N. Venisse, S. Marchand, M. Bachelet, W. Couet, A simple and sensitive
781 liquid chromatography-tandem mass spectrometry assay for the quantification of ertapenem
782 in microdialysate, *Journal of chromatography. B, Analytical technologies in the biomedical*
783 *and life sciences* 862(1-2) (2008) 242-5.

784 [47] G. la Marca, E. Giocaliere, F. Villanelli, S. Malvagia, S. Funghini, D. Ombrone, L.
785 Filippi, M. De Gaudio, M. De Martino, L. Galli, Development of an UPLC-MS/MS method
786 for the determination of antibiotic ertapenem on dried blood spots, *Journal of pharmaceutical*
787 *and biomedical analysis* 61 (2012) 108-13.

788 [48] T. Ohmori, A. Suzuki, T. Niwa, H. Ushikoshi, K. Shirai, S. Yoshida, S. Ogura, Y. Itoh,
789 Simultaneous determination of eight beta-lactam antibiotics in human serum by liquid
790 chromatography-tandem mass spectrometry, *J Chromatogr B* 879(15-16) (2011) 1038-42.

791 [49] M.J. Schneider, S.J. Lehotay, A.R. Lightfield, Evaluation of a multi-class, multi-residue
792 liquid chromatography-tandem mass spectrometry method for analysis of 120 veterinary
793 drugs in bovine kidney, *Drug testing and analysis* 4 Suppl 1 (2012) 91-102.

794 [50] L. Geis-Asteggianti, S.J. Lehotay, A.R. Lightfield, T. Dutko, C. Ng, L. Bluhm,
795 Ruggedness testing and validation of a practical analytical method for >100 veterinary drug
796 residues in bovine muscle by ultrahigh performance liquid chromatography–tandem mass
797 spectrometry, *Journal of Chromatography A* 1258 (2012) 43-54.

798 [51] B.J. Berendsen, I.J. Elbers, A.A. Stolker, Determination of the stability of antibiotics in
799 matrix and reference solutions using a straightforward procedure applying mass spectrometric
800 detection, *Food Additives & Contaminants: Part A* 28(12) (2011) 1657-66.

801 [52] L. Okerman, J. Van Hende, L. De Zutter, Stability of frozen stock solutions of beta-
802 lactam antibiotics, cephalosporins, tetracyclines and quinolones used in antibiotic residue
803 screening and antibiotic susceptibility testing, *Analytica Chimica Acta* 586(1-2) (2007) 284-
804 288.

805 [53] K.L. Tyczkowska, R.D. Voyksner, A.L. Aronson, Solvent degradation of cloxacillin in
806 vitro, *Journal of Chromatography A* 594(1) (1992) 195-201.

807 [54] S. Pellicciotti, S. Moretti, G. Saluti, R. Galarini, G. Dusi, Stability of antibiotics in
808 solution: a critical issue during the development of a multi-class method by LC-HRMS
809 (Orbitrap), 7th International Symposium On Hormone and Veterinary Drug Residues
810 Analysis, Ghent (Belgium), 2–5 June 2014.

811 [55] W.A. Moats, K.L. Anderson, J.E. Rushing, S. Buckley, Conversion of cephapirin to
812 deacetylcephapirin in milk and tissues of treated animals, *Journal of Agricultural and Food*
813 *Chemistry* 48(2) (2000) 498-502.

814 [56] D.N. Heller, D.A. Kaplan, N.G. Rummel, J. von Bredow, Identification of cephalirin
815 metabolites and degradants in bovine milk by electrospray ionization-ion trap tandem mass
816 spectrometry, J. Agric. Food Chem. 48(12) (2000) 6030-6035.

817

Table 1. UHPLC-MS/MS conditions for β -lactam antibiotics.

Analyte	Measured ion	RT (min)	Precursor ion (m/z)	Product ions (m/z)	Dwell time (s)	CV (V)	CE (eV)	MRM ^b
Imipenem	[M+H] ⁺	2.28	300.1	142.0 ^a /126.1	0.194	22	28/17	1
Biapenem	[M+H] ⁺	2.81	351.2	110.2 ^a /170.1	0.180	25	19/20	2
Amoxicillin	[M+H] ⁺	3.21	366.1	114.0 ^a /349.1	0.035	18	20/9	3
Amoxicillin-D ₄	[M+H] ⁺	3.21	370.6	354.0	0.020	16	9	3
Cefadroxil	[M+H] ⁺	3.24	364.1	208.1 ^a /158.0	0.030	15	9/9	3
Cefadroxil-D ₄	[M+H] ⁺	3.23	368.6	115.0	0.012	16	20	3
DAC	[M+H] ⁺	3.28	382.1	152.1 ^a /226.1	0.030	29	24/23	3
DAC-D ₆	[M+H] ⁺	3.27	387.9	115.1	0.015	31	45	3
Doripenem	[M+H] ⁺	3.35	421.1	342.2 ^a /274.2	0.021	20	13/16	3
Meropenem	[M+H] ⁺	3.50	384.1	141.0 ^a /114.1	0.103	23	15/23	4
Ampicillin	[M+H] ⁺	3.75	350.1	106.1 ^a /192.1	0.033	22	21/16	5
Ampicillin-D ₅	[M+H] ⁺	3.74	355.0	197.1	0.025	18	16	5
Cephalexin	[M+H] ⁺	3.74	348.1	158.1 ^a /174.1	0.036	20	9/17	5
Cephalexin-D ₅	[M+H] ⁺	3.73	353.0	179.1	0.025	17	14	5
Cephapirin	[M+H] ⁺	3.76	424.0	292.1 ^a /152.1	0.040	24	15/24	5
Cefquinome	[M+H] ⁺	4.05	529.0	134.2 ^a /167.1	0.074	21	13/23	6
Ertapenem	[M+H] ⁺	4.19	476.1	432.1 ^a /233.2	0.016	20	9/14	6
Mecillinam	[M+H] ⁺	4.19	326.1	167.2 ^a /139.2	0.022	34	15/29	6
DCCD	[M+H] ⁺	4.26	548.7	183.1 ^a /241.2	0.032	36	29/20	6
DCCD-D ₃	[M+H] ⁺	4.26	551.9	366.0	0.032	36	21	6
Cefalonium	[M+H] ⁺	4.24	458.9	152.0 ^a /337.0	0.006	18	18/10	6
Cefacetile	[M+NH ₄] ⁺	4.92	357.0	280.1 ^a /252.1	0.056	16	9/14	7
Cefotaxime	[M+H] ⁺	4.92	456.1	396.0 ^a /167.1	0.080	23	10/19	7
Faropenem	[M+Na] ⁺	5.04	308.0	178.0 ^a /222.1	0.047	26	18/12	8
Cefazolin	[M+H] ⁺	5.13	454.9	156.0 ^a /295.1	0.047	17	15/15	8
Cefazolin- ¹³ C ₂ ¹⁵ N	[M+H] ⁺	5.13	457.8	298.1	0.047	15	16	8
Cefuroxime	[M+NH ₄] ⁺	5.39	442.0	364.1 ^a /336.1	0.129	16	14/14	9
Cefoperazone	[M+H] ⁺	5.57	646.0	143.2 ^a /530.0	0.022	25	33/12	10
DCD	[M+2H] ²⁺	5.71	429.2	183.1 ^a /397.0	0.026	23	21/13	10
Ticaracillin	[M+H] ⁺	5.82	385.1	160.1 ^a /114.1	0.039	25	14/36	10
Methicillin	[M+H] ⁺	5.98	381.1	165.2 ^a /222.2	0.038	27	23/18	10
Piperacillin	[M+H] ⁺	6.02	518.0	143.2 ^a /160.1	0.020	28	17/10	10
Penicillin G	[M+H] ⁺	6.23	335.1	160.1 ^a /176.1	0.135	21	12/12	11
Penicillin G-D ₇	[M+H] ⁺	6.21	342.0	183.1	0.080	18	11	11
Penicillin V	[M+H] ⁺	6.54	351.1	160.1 ^a /114.1	0.139	21	13/33	12
Penicillin V-D ₅	[M+H] ⁺	6.52	355.9	160.0	0.050	20	15	12
Oxacillin	[M+H] ⁺	6.76	402.1	160.1 ^a /243.2	0.150	20	16/14	13
Cloxacillin	[M+H] ⁺	7.04	436.0	160.1 ^a /277.1	0.115	22	12/12	14
Nafcillin	[M+H] ⁺	7.16	415.1	199.2 ^a /171.2	0.081	20	7/28	14
Nafcillin-D ₅	[M+H] ⁺	7.15	420.0	261.1	0.050	19	17	14
Dicloxacillin	[M+H] ⁺	7.41	470.0	160.2 ^a /114.2	0.151	22	14/34	15

a = Quantitation ion; RT = Retention Time; CV = Cone Voltage; CE = Collision Energy.

b = MRM Windows: 1 = (1.99 – 2.98 min); 2 = (2.60 – 3.13 min); 3 = (2.98 – 3.50 min); 4 = (3.41 – 3.62 min); 5 = (3.60 – 3.88 min); 6 = (3.91 – 4.42 min); 7 = (4.75 – 5.10 min); 8 = (4.93 – 5.26 min); 9 = (5.23 – 5.57 min); 10 = (5.35 – 6.27 min); 11 = (6.05 – 6.36 min); 12 = (6.36 – 6.69 min); 13 = (6.62 – 6.99 min); 14 = (6.82 – 7.35 min); 15 = (7.23 – 7.58 min).

Table 2. Calibration ranges, LODs, LOQs and matrix effect data. Positive values indicate ion suppression, while negative values indicate matrix enhancement.

Analyte	Calibration range ($\mu\text{g kg}^{-1}$)	R^2	LOD ($\mu\text{g kg}^{-1}$)	LOQ ($\mu\text{g kg}^{-1}$)	ME% ^a	RSD(%) (no IS)	RSD(%) (IS)
Amoxicillin	12.5-100	0.994	0.50	1.7	75	17.2	4.6
Ampicillin	12.5-100	0.994	0.050	0.17	34.1	12.2	3.2
Cloxacillin	75-600	0.998	0.15	0.50	13.8	4.7	4.5
Dicloxacillin	75-600	0.998	0.10	0.30	10.6	6.7	5.5
Mecillinam	2.5-20	0.996	0.030	0.10	19.7	27.0	26.2
Methicillin	2.5-20	0.994	0.030	0.10	-73	19.6	No IS used
Nafcillin	75-600	0.998	0.030	0.10	9.5	4.2	1.6
Oxacillin	75-600	0.992	0.010	0.035	17.0	5.4	4.9
Penicillin G	12.5-100	0.993	0.040	0.12	19.2	7.2	4.1
Penicillin V	6.25-50	0.998	0.020	0.070	14.3	3.6	3.7
Piperacillin	2.5-20	0.997	0.040	0.12	5.8	9.3	No IS used
Ticarcillin	5-40	0.993	0.30	1.0	21.7	10.0	No IS used
Cefacetrile	5-40	0.994	0.30	1.0	44.6	22.7	12.5
Cefadroxil	12.5-100	0.988	1.3	4.2	69	22.0	11.0
Cephalexin	50-400	0.998	0.50	1.7	23.7	10.9	3.2
Cefalonium	5-40	0.989	0.30	1.0	42.0	20.2	13.9
Cefazolin	2.5-20	0.997	0.10	0.30	35.8	13.8	3.7
Cefoperazone	5-40	0.994	0.30	1.0	-15.2	9.2	No IS used
Cefotaxime	5-40	0.994	0.30	1.0	43.3	24.6	14.1
Cefquinome	12.5-100	0.989	0.070	0.22	-34.1	28.9	No IS used
Cefuroxime	5-40	0.995	0.30	1.0	6.8	7.8	7.7
DAC	12.5-100	0.995	0.60	1.9	70	23.6	5.8
DCCD	125-1000	0.993	4.8	16.0	44.1	14.3	4.6
DCD	125-1000	0.990	0.20	0.70	52	10.1	No IS used
Biapenem	5-40	0.992	0.30	1.0	-34.8	19.5	No IS used
Doripenem	12.5-100	0.991	1.0	3.0	42.9	25.0	21.6
Ertapenem	50-400	0.994	1.5	5.0	30.5	8.6	8.7
Imipenem	12.5-100	0.994	0.50	1.7	-45.8	15.2	No IS used
Meropenem	5-40	0.997	0.30	1.0	61	12.3	No IS used
Faropenem	12.5-100	0.995	0.50	1.7	48.1	35.8	25.0

^a = Results are based on the analysis of 27 different bovine muscle samples; IS = Internal Standard.

Table 3. Validation results for bovine muscle.

Analyte	MRL/TL ($\mu\text{g kg}^{-1}$)	Trueness% (WLR%)			Trueness% (WLR%)			CC α ($\mu\text{g kg}^{-1}$)	CC β ($\mu\text{g kg}^{-1}$)	OAR% (RSD%, n = 3)	
		0.5 MRL/TL	MRL/TL	1.5 MRL/TL	0.5 MRL/TL	MRL/TL	1.5 MRL/TL			2LCL	0.875HCL
Amoxicillin	50	102 (5.9)	101 (4.5)	103 (5.2)	101 (6.4)	101 (5.4)	101 (5.8)	55	60	69 (2.4)	70 (5.5)
Ampicillin	50	101 (2.4)	99 (2.4)	101 (3.3)	100 (4.3)	102 (3.6)	100 (4.0)	54	58	78 (7.9)	76 (4.8)
Cloxacillin	300	103 (2.2)	101 (2.2)	101 (1.7)	102 (6.0)	99 (5.2)	95 (7.5)	350	399	81 (5.6)	88 (3.6)
Dicloxacillin	300	104 (2.4)	103 (2.3)	104 (5.0)	102 (6.2)	99 (3.7)	97 (5.6)	321	351	83 (2.3)	87 (3.4)
Mecillinam	10	96 (2.5)	97 (3.3)	100 (4.0)	81 (29.0)	69 (29.9)	91 (21.6)	13.1	19.6	85 (1.1)	83 (2.5)
Methicillin	10	108 (10.2)	104 (9.5)	98 (6.7)	91 (13.2)	90 (14.8)	91 (13.0)	12.2	15.8	85 (5.0)	84 (4.4)
Nafcillin	300	101 (2.2)	101 (1.8)	100 (1.5)	104 (6.2)	101 (3.2)	99 (4.2)	326	350	84 (2.0)	88 (4.2)
Oxacillin	300	103 (2.7)	101 (3.2)	101 (3.2)	105 (6.4)	103 (3.6)	101 (4.3)	322	347	78 (12.6)	89 (5.8)
Penicillin G	50	106 (2.8)	105 (2.8)	105 (2.4)	100 (6.7)	104 (12.6)	103 (14.9)	62	84	79 (12.5)	85 (3.1)
Penicillin V	25	101 (2.4)	100 (2.3)	102 (2.8)	98 (4.8)	98 (2.3)	100 (2.8)	26.0	27.3	87 (1.5)	86 (3.4)
Piperacillin	10	104 (4.9)	104 (5.8)	105 (5.3)	102 (7.4)	102 (6.1)	102 (5.3)	11.5	12.6	89 (4.0)	85 (4.6)
Ticarillin	20	99 (6.4)	99 (8.6)	100 (6.7)	99 (10.0)	105 (10.3)	95 (8.9)	28.5	34.1	61 (5.3)	62 (6.6)
Cefacetrile	20	101 (4.3)	100 (4.6)	99 (8.2)	108 (9.6)	104 (5.5)	101 (8.9)	22.9	27.1	80 (2.7)	79 (3.4)
Cefadroxil	50	100 (5.7)	98 (5.4)	101 (6.3)	101 (10.9)	98 (5.6)	97 (6.7)	54	60	58 (2.1)	61 (8.1)
Cephalexin	200	98 (1.9)	98 (1.9)	100 (3.5)	101 (3.5)	100 (2.0)	99 (3.6)	202	213	69 (4.8)	67 (3.9)
Cefalonium	20	100 (8.9)	100 (5.6)	100 (6.0)	99 (12.6)	98 (8.0)	101 (8.5)	21.9	25.0	79 (13.0)	81 (2.4)
Cefazolin	10	101 (5.0)	101 (4.5)	100 (5.2)	98 (6.8)	101 (5.6)	101 (8.3)	10.6	11.9	81 (4.6)	83 (4.2)
Cefoperazone	20	104 (5.5)	101 (7.0)	102 (6.3)	101 (6.7)	100 (8.5)	101 (8.0)	23.7	27.4	86 (3.7)	86 (6.7)
Cefotaxime	20	102 (6.4)	104 (5.2)	101 (6.5)	101 (6.6)	104 (6.1)	107 (8.8)	21.9	25.0	64 (5.0)	67 (7.0)
Cefquinome	50	103 (7.0)	101 (11.5)	100 (6.9)	95 (10.5)	96 (12.7)	91 (7.1)	68	79	87 (3.5)	86 (5.3)
Cefuroxime	20	104 (7.2)	102 (6.7)	101 (7.2)	100 (7.2)	105 (6.7)	104 (10.7)	21.1	24.7	87 (1.4)	84 (3.2)
DAC	50	102 (4.9)	101 (4.4)	101 (3.7)	99 (6.8)	97 (5.2)	99 (4.1)	53	57	77 (5.0)	78 (2.7)
DCCD	500 ^a	104 (4.8)	101 (5.9)	104 (4.7)	101 (4.8)	99 (5.9)	99 (4.7)	536	584	40.3 (7.8)	45.1 (2.9)
DCD	500 ^a	108 (5.2)	102 (7.0)	103 (8.8)	99 (9.7)	106 (9.4)	98 (8.8)	677	802	53 (4.7)	56 (7.9)
Biapenem	20	106 (5.6)	103 (2.9)	103 (5.8)	91 (14.2)	84 (13.9)	92 (8.2)	23.0	27.3	71 (2.2)	71 (3.6)
Doripenem	50	106 (4.8)	103 (7.1)	100 (4.7)	96 (10.9)	94 (8.8)	96 (12.3)	57	71	62 (3.9)	65 (2.8)
Ertapenem	200	101 (4.7)	97 (4.3)	102 (8.0)	96 (9.0)	95 (7.3)	100 (8.4)	213	242	44.2 (6.4)	43.3 (6.8)
Imipenem	50	106 (6.0)	100 (8.4)	86 (10.9)	96 (10.9)	102 (10.2)	86 (13.9)	69	90	47.0 (8.3)	51 (4.9)



Table 3. *(Continued)*

Analyte	MRL/TL ($\mu\text{g kg}^{-1}$)	Trueness% (WLR%)			Trueness% (WLR%)			CC α ($\mu\text{g kg}^{-1}$)	CC β ($\mu\text{g kg}^{-1}$)	OAR% (RSD%, n = 3)	
		0.5 MRL/TL	MRL/TL	1.5 MRL/TL	0.5 MRL/TL	MRL/TL	1.5 MRL/TL			2LCL	0.875HCL
Meropenem	20	104 (3.6)	100 (6.0)	99 (3.6)	115 (21.5)	125 (23.4)	122 (22.4)	28.4	44.1	66 (0.46)	66 (8.4)
Faropenem	50	97 (17.3)	96 (14.8)	97 (11.7)	128 (17.8)	143 (14.9)	119 (13.5)	77	100	78 (0.10)	79 (3.3)

TL = Target Level; OAR = Overall Absolute Recovery; LCL = Lowest Calibration Level; HCL = Highest Calibration Level.

a = Validation levels were 250, 500 and 750 $\mu\text{g kg}^{-1}$ (MRL established by (EU) 37/2010 is 1000 $\mu\text{g kg}^{-1}$ for the sum of all residues retaining the β -lactam structure, expressed as desfuroylceftiofur).



Penicillins




Mecillinam

<i>Penicillin</i>	<i>R</i>	<i>Penicillin</i>	<i>R</i>
Amoxicillin	-C ₇ H ₈ NO	Oxacillin	-C ₁₀ H ₈ NO
Ampicillin	-C ₇ H ₈ N	Penicillin G	-C ₇ H ₇
Cloxacillin	-C ₁₀ H ₇ ClNO	Penicillin V	-C ₇ H ₇ O
Dicloxacillin	-C ₁₀ H ₆ Cl ₂ NO	Piperacillin	-C ₁₄ H ₁₆ N ₃ O ₃
Methicillin	-C ₈ H ₉ O ₂	Ticarcillin	-C ₆ H ₅ O ₂ S
Nafcillin	-C ₁₂ H ₁₁ O		

Cephalosporins





Cefalonium
and
Cefquinome




<i>Cephalosporin</i>	<i>R1</i>	<i>R2</i>	<i>Cephalosporin</i>	<i>R1</i>	<i>R2</i>
Cefacetrile	-C ₃ H ₅ O ₂	-C ₂ H ₂ N	Cefotaxime	-C ₃ H ₅ O ₂	-C ₅ H ₆ N ₃ OS
Cefadroxil	-CH ₃	-C ₇ H ₈ NO	Cefquinome	-C ₁₀ H ₁₃ N	-C ₅ H ₆ N ₃ OS
Cephalexin	-CH ₃	-C ₇ H ₈ N	Cefuroxime	-C ₂ H ₄ NO ₂	-C ₆ H ₆ NO ₂
Cefalonium	-C ₇ H ₈ N ₂ O	-C ₅ H ₅ S	DAC	-CH ₃ O	-C ₆ H ₆ NS
Cefazolin	-C ₄ H ₅ N ₂ S ₂	-C ₂ H ₃ N ₄	DCD	-CH ₂ S	-C ₅ H ₆ N ₃ OS
Cefoperazone	-C ₃ H ₅ N ₄ S	-C ₁₄ H ₁₆ N ₃ O ₄	DCCD	-C ₄ H ₈ NO ₂ S ₂	-C ₅ H ₆ N ₃ OS

Carbapenems and faropenem



Biapenem



Faropenem

<i>Carbapenem</i>	<i>R1</i>	<i>R2</i>
Biapenem	-CH ₃	-C ₅ H ₇ N ₃ S
Doripenem	-CH ₃	-C ₅ H ₁₂ N ₃ O ₂ S ₂
Ertapenem	-CH ₃	-C ₁₂ H ₁₃ N ₂ O ₃ S
Imipenem	-H	-C ₃ H ₇ N ₂ S
Meropenem	-CH ₃	-C ₇ H ₁₃ N ₂ OS

Fig. 1. Chemical structures of the β -lactams included in the proposed UHPLC-MS/MS method.





Fig. 2. Overlay of LC-MS/MS chromatograms for the β -lactams (A-C) and internal standards (D-E) included in the method in a blank bovine muscle sample fortified at half the MRL/TL.





Fig. 2. (Continued)




Fig. 2. (*Continued*)




Fig. 3. Average recoveries (and standard deviations, shown by error bars) obtained using NH₂, PSA, Z-Sep and C₈ sorbents or hexane liquid-liquid partitioning clean-up (n = 3).




Fig. 4. Average recoveries (and standard deviations, shown by error bars) obtained using C₁₈, Z-Sep+ and Z-Sep/C₁₈ sorbents (n = 3).





Fig. 5. MS spectra obtained from a cloxacillin standard solution prepared in $\text{H}_2\text{O}:\text{MeCN}$ (A) and a cloxacillin standard solution prepared in MeOH (B), after 11 days of storage at $+4^\circ\text{C}$.




Fig. 6. Extracted chromatograms for the masses 436.0 m/z (A) and 468.0 m/z (B) obtained after column injections of 10 $\mu\text{g mL}^{-1}$ cloxacillin standard solutions prepared in H₂O:MeCN or MeOH and stored for 11 days at different temperatures. RT = Room Temperature.




Fig. 7. Flow diagram of published methods (Becker *et al.* (2004) [41]; Berendsen *et al.* (2013) [2]) sample preparation procedures and comparison with the proposed method for the extraction of β -lactam residues from muscle tissue samples, showing steps and time required. The estimated time for each step is based on the extraction of one sample only.